

# MAST URI<sup>®</sup> CONNECT

The Connected Solution for Urine Microbiology



Direct Susceptibility  
Testing - Results in  
18 hours



Phenotypic Results  
Analysis - ID, AST  
and AMR



End-to-End Sample  
Traceability



Integrated EUCAST-  
aligned Reading  
Rules

Urinary Tract Infections (UTIs) are one of the most common conditions presenting in primary healthcare. Traditional methods of urine culturing are labour intensive and slow. Laboratories are often left facing rising sample numbers with the expectation of faster turnaround times, often without additional resources. The **MAST URI® CONNECT** provides the connected solution to Urine Microbiology.



**MAST URI®**  
**WELL CONNECT**



**MAST URI®**  
**PREP CONNECT**



**MAST URI®**  
**AUTOPREP CONNECT**

## Sample Preparation

**MAST URI® CONNECT** provides laboratories with the opportunity to select from three sample dispensing options. Consistent with all options is the use of a barcoded template, a Urine Specimen Plate (USP). This bridges the gap in sample traceability, connecting the barcoded patient samples are dispensed in a USP, using one of the 3 options.

### **MAST URI® WELL CONNECT**

A semi-automated dispensing aid that simplifies workflow, initiates end-to-end traceability and aids in accurate sample dispensing.

### **MAST URI® PREP CONNECT**

Utilising the **MAST URI® WELL CONNECT**, this option allows for a user to perform sample preparation while simultaneously freeing up the **MAST URI® PLUS CONNECT** for plate reading and result validation.

### **MAST URI® AUTOPREP CONNECT**

A fully automated, walkaway dispensing solution that frees up valuable laboratory team members, making this ideal for high throughput requirements.

## Sample Inoculation

Following the completion of the dispensing process, the Urine Specimen Plate (USP) is then used to inoculate bespoke plate sets using multipoint inoculation. **MAST URI® DOT CONNECT** utilises barcode readers for both the USP template and **MAST URI® Plates** to facilitate end-to-end traceability, and tracks pin sterilisation through audited disinfection cycles.



**MAST URI®**  
**DOT CONNECT**



**MAST URI®**  
**PLUS CONNECT**


## Sample Capture

After incubation, the **MAST URI® Plates** are imaged on the **MAST URI® PLUS CONNECT**, with images and sample information being imported into the database for validation. The **MAST URI® PLUS CONNECT** saves precious bench space by incorporating a built-in PC, boasts a high-resolution camera for sharper image capture and possesses refined software for seamless sample analysis.

## Sample Validation

The overhauled **MAST URI® CONNECT** software has been built from the ground up, focusing on greater accessibility, streamlining workflow and having a strong focus on traceability. It encompasses EUCAST expert reading rules, assisting in the exporting of clinically relevant results. The combination of sharper images and an all-encompassing sample analysis screen, ensures an efficient validation process presenting a complete collection of sensitivity and identification.

Once validated, results are exported to LIMS, as the **MAST URI® CONNECT** can be fully interfaced.



# Next Generation Urine Analysis Get **CONNECT**ed

Product development based on real customer feedback.

*“Having to remove the pin head after every batch could increase contamination risk”*

**MAST URI® CONNECT** has integrated a mandatory sterilisation step between inoculation of plate sets, which aids in preventing cross-contamination

*“Increasing sample traceability would improve auditability”*

End-to-end traceability is a fundamental aspect of the **MAST URI® CONNECT** with integrated barcode scanners enabling seamless information transfer between the whole suite of instruments

*“It would be great to be able to export results for urgent samples”*

With advanced software, **MAST URI® CONNECT** can now export single samples providing increased flexibility and allowing urgent results to be exported faster

**MAST URI®**  
**CONNECT**

The enhancements to the user interface on the **MAST URI® CONNECT** allows for manual entry of typed notes at both sample dispensing and result validation

*“Manually writing notes can be time consuming and lead to errors”*

Integrated fail-safes, including an incubation warning signal, decreases likelihood of user error and improves workflow consistency across the whole system


*“More guidance through the workflow would be helpful for new team members”*

**MAST URI® CONNECT** boasts a higher resolution camera providing high quality images allowing for well images during result validation


*“Having better images would make result validation easier”*

# The MAST URI<sup>®</sup> CONNECT reduces staff processing time by up to 70%<sup>1</sup>

**MAST URI<sup>®</sup> CONNECT process > 18 - 24 hours**



**Traditional process > 36 - 48 hours**




**MAST URI<sup>®</sup> CONNECT** is a sustainable alternative to the traditional process with an **85%** reduction in plastic and **96%** reduction in agar waste.<sup>2</sup>


Our **MAST URI<sup>®</sup> Plates** are produced according to current EUCAST breakpoint guidelines.

The new **MAST URI<sup>®</sup> CONNECT** builds on its renowned predecessor the **MAST URI<sup>®</sup> System**.

By using fewer consumables **MAST URI<sup>®</sup> CONNECT** streamlines urine AST and saves on storage space.



Samples incubated




AST profile obtained



Results validated

**References:**

1. Pathology in Practice. Streamlining urine microbiology: the Northampton experience. Vol 13(3) 2012.
2. Mast Group Ltd. Study conducted comparing 14 MPM plate set to 96 Petri dishes. 2023.



# Introducing **MAST URI<sup>®</sup> AUTOPREP CONNECT**



### Specimen Tube Compatibility

A range of tubes commonly used for urine sample collection can be processed on the device.



### Easy Setup

Internal deck layout is designed for rapid loading and easy cleaning for a faster workflow.



### Barcode Recognition

Integrated barcode readers allow the device to maintain full end-to-end traceability for tracking of samples.



### Automatic Decapping


Decaps all samples reducing the likelihood of cross-contamination from manual handling and reduces risk of RSI.



### Vortexing and Pipetting


Automatic vortexing and liquid level detection of each sample occurs before sample aspiration and dispensing.

## New **MAST URI<sup>®</sup> CONNECT** Interface



### Improved User Interface

Upgraded software technology and a new LED touch screen monitor provides customers with improved accessibility.



### Simplified Workflow Process

Simultaneously display ID, sensitivity data, patient demographics, EUCAST reading rules and manual typed notes.



### High quality well images

Higher resolution images with the new zoom function greatly enhances efficiency for validation, further reducing reporting times.




# MAST URI® CONNECT

## Product Range

### Product Specifications

MAST URI® CONNECT full product range.

#### MAST URI® AUTOPREP CONNECT



#### MAST URI® WELL CONNECT



#### MAST URI® PREP CONNECT



<b>Weight</b>	200 kg	1.2 kg	Tablet: 0.80 kg USB Hub: 0.52 kg Well Connect: 1.2 kg
<b>Dimensions</b>	1450W × 650D × 1088H mm	300W × 150D × 80H mm	Tablet: 292W × 8.5D × 201H mm USB Hub: 130W × 30D × 70H mm Well: 300W × 150D × 80H mm
<b>Power requirements</b>	240 VAC, 50/60 Hz, 500 W	5 V, 500 mA, USB powered	240 VAC, 50/60 Hz, 199 W
<b>Connectivity</b>	Ethernet, USB	–	Ethernet, USB
<b>Product code</b>	URIAUTOPREPCONNECT*	URIWELLCONNECT	URIPREPCONNECT

#### MAST URI® DOT CONNECT



#### MAST URI® PLUS CONNECT



#### MAST URI® PLUS CONNECT *Hardware*



<b>Weight</b>	11.5 kg	21.5 kg	3.9 kg
<b>Connectivity</b>	380W × 415D × 290H mm	522W × 522D × 671H mm	565W × 175D × 433H mm
<b>Power requirements</b>	100-240 VAC, 50/60 Hz, 60 W	100-240 VAC, 50/60 Hz, 220 W	90-264 VAC, 50-60 Hz, 25 W
<b>Connectivity</b>	Ethernet	USB, Ethernet, Display Port	USB C, Ethernet
<b>Product code</b>	URIDOTCONNECT	URIPLUSCONNECT*	URIPLUSCONNECT

\* These devices will be fitted with an uninterruptible power supply (UPS)

# MAST URI® Plates

**MAST URI® Plates:** Antibiotic and Identification plates; current availability (subject to change)

Antibiotic plates aligned with EUCAST Clinical Breakpoints V15.0


Code	Plate Type
MPM-AK8	Amikacin 8mg/L
MPM-AMX8	Amoxicillin 8mg/L <sup>1</sup>
MPM-AMXC8	Amoxicillin-Clavulanate 8mg/L <sup>1</sup>
MPM-AMXC32	Amoxicillin-Clavulanate 32mg/L <sup>1</sup>
MPM-AMP8	Ampicillin 8mg/L <sup>1</sup>
MPM-CFX16	Cefalexin 16mg/L
MPM-CXM8	Cefuroxime 8mg/L
MPM-CPD1	Cefpodoxime 1mg/L
MPM-CIP0.25	Ciprofloxacin 0.25mg/L
MPM-CIP0.5	Ciprofloxacin 0.5mg/L
MPM-CRO1	Ceftriaxone 1mg/L

Code	Plate Type
MPM-FOS8	Fosfomycin 8mg/L
MPM-GM2	Gentamicin 2mg/L
MPM-MEC8	Mecillinam 8mg/L
MPM-NI64	Nitrofurantoin 64mg/L
MPM-PTZ8	Piperacillin-Tazobactam 8mg/L
MPM-TEM16	Temocillin 16mg/L
MPM-TM4	Trimethoprim 4mg/L
MPM-TS2/38	Trimethoprim 2mg/L + Sulfamethoxazole 38mg/L
MPM-VA4	Vancomycin 4mg/L

1. For information on how to implement aminopenicillin breakpoints, see [www.eucast.org/eucastguidancedocuments/](http://www.eucast.org/eucastguidancedocuments/)

## Identification Plate Set

Code	Description
MPM-CON	Control
MPM-CPOS	Chromogenic Gram-positive
MPM-CNEG	Chromogenic Gram-negative
MPM-CBGA	Chromogenic Beta Glucuronidase
MPM-IDM51	Aesculin
MPM-IDM52	Phenylpyruvic Acid
MPM-IDM62	Chromogenic Group B Strep



MPM-CON

## Detection of Antimicrobial Activity

Code	Description
MPM-IDM54	Antimicrobial Activity Indicator (AAI) plate

## Plates for Antimicrobial Resistance Detection

Code	Description
MPM-CRE	Carbapenemase screening plate
MPM-CPD4	Detection of ESBLs. Eliminates bacteria lying on EUCAST breakpoint (1mg/L)

### Notes:

Direct culture and sensitivity testing of samples from patients on antibiotics may give misleading results due to inhibition of pathogens. To facilitate detection of antimicrobial activity in urine samples the Antimicrobial Activity Indicator plate (MPM-IDM54) may be included in the Plate Set if required.

**These plates are currently available with a four-week lead-time.**

### United Kingdom

Mast Group Ltd.  
Mast House, Derby Road, Bootle  
Merseyside L20 1EA

+44 (0)151 933 7277  
sales@mastgrp.com

### Germany

Mast Diagnostica GmbH  
Feldstraße 20  
23858 Reinfeld

+49 (0)4533 2007 0  
mast@mast-diagnostica.de

### France

Mast Diagnostic  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1

+33 (0)322 80 80 67  
info@mast-diagnostic.fr

### South Africa

Davies Diagnostics  
145 Oak Avenue, Ferndale  
Randburg, Gauteng, 2194

+27 11 777 7600  
info@daviesdiag.co.za